


REV	ECN	Revision History	Date	Editor
A05	D00021277	Software updates to V3.0.7.2.	2015-07-07	Chen Ying
A04		Software updates to V3.0.7.1.	2015-03-23	Lei Wanyi
A03	D00018915	Change address of manufacturer	2015-03-19	Lei Wanyi
A02	D00016369	Software updates to V3.0.6.1.	2014-10-13	Lei Wanyi
A01	D00008056	First release	2013-07-10	Xu Yinghua

Title	S9 Technical Specifications
--------------	------------------------------------

This document contains intellectual property information that is proprietary to SonoScape Medical Corp. and is protected by law. Neither the document nor the information contained therein should be used or reproduced in whole or partially, without prior written agreement consent of SonoScape Medical Corp.

 SonoScape Medical Corp.	Document Number 901-00724		Distribution Number
	Version A05	Effective Date 2015-07-07	Page Page 1 of 32

Specifications for S9 Portable Digital Color Doppler Ultrasound System



SonoScape

THE PIONEER OF DOPPLER ULTRASOUND IN CHINA

Overview

The S9 Color Doppler Ultrasound System adopts the advanced ultrasonic Doppler technologies, including the Full Digital Super-wide Band Beam Former, Digital Dynamic Focusing, Variable Aperture and Dynamic Tracing, Wide Band Dynamic Range, Multi-Beam Parallel Processing, etc. The ultrasound diagnostic software system, ultrasound system imaging, multi-languages operation interfaces and touch screen with human-computer interaction technology can be customized easily in accordance with the design of human engineering. Users can operate the system with the minimum requirement of training or guidance. This system has been designed to comply with applicable international standards and regulations, ensuring the safety and availability of this product.

This system is based on the computer technology and Linux operation system, which make the system more flexible and stable. System maintenance and function update can be completed by software updating, which promotes product value and keeps the technological advancement.

Advanced Technologies

- New generation digital front-end technology
- Multi-beam processing technology
- Spatial compound imaging
- μ Scan image processing technology
- Tissue Harmonic Imaging
- Invert harmonic imaging
- High Pulse Repetition Frequency
- Panoramic imaging
- 4D imaging
- Graphic diagnosis icon

- Touch screen with human-computer interaction technology
- Elastography imaging
- Contrast imaging

Standard Configurations

- Touch screen
- ECG function module
- μ Scan function
- B mode: Five variable frequencies
- LGC (two levels)
- Tissue acoustic characteristics index
- THI mode
- PIH mode
- Multi-beam
- Image rotation
- Compound imaging
- Trapezoidal imaging
- Color mode
- PDI mode
- DPDI mode
- PW mode
- HPRF support
- Simult mode
- FreeHand 3D
- Steer M mode
- Color M function
- TDI function
- CW mode
- Stress Echo function
- 2D panoramic imaging
- Color flow panoramic imaging
- Biopsy enhance
- Basic measurement package

- Obstetrics measurement package
- Gynecology measurement package
- Cardiology measurement package
- Abdomen measurement package
- Vascular measurement package
- Urology measurement package
- Small parts measurement package
- Pediatrics measurement package
- TEI index
- PW auto trace
- IMT measurement
- DICOM transmission
- DICOM Worklist
- DICOM MPPS
- DICOM C-store
- DICOM Q/R

Optional Functions

- 4D imaging
- Biplane dual live
- Elastography imaging
- Contrast imaging
- DVD

Optional Accessories

- Biopsy bracket
- Color ink-jet printer
- B/W video printer
- Color video printer
- Probe cable hanger
- Foot switch
- External DVD burner
- Trolley

Scan Methods

- Electronic curved sector scan
- Electronic linear array scan
- Electronic phased array sector scan
- Mechanical sector scan (convex probe)

Applications

- Abdomen
- Vascular
- Cardiology
- Gynecological/Obstetrics
- Urology
- Musculoskeletal
- Interventional ultrasound
- Small parts
- Anesthesiology
- Pediatrics
- Orthopedic
- Cephalic

Imaging Modes

- B mode
- M mode
- THI mode
- CFM mode
- PDI mode
- TDI mode
- PW mode
- CW mode
- HPRF
- 3D/4D mode
- Color M mode
- Steer M mode

- Elastography imaging
- Contrast imaging

Display Formats

- Dual B
- Quad B
- B + Color
- B + Color dual live
- B + PDI
- B + TDI
- B + M
- B + Color + M
- B + PDI + M
- B + TDI + M
- B + PW
- B + Color + PW
- B + PDI + PW
- B + TDI + PW
- B + CW
- B + Color + CW
- B + PDI + CW
- Panoramic imaging
- Trapezoidal imaging

System Settings Menu

- General
 - General
 - Hospital Name
 - Language
 - ⊕ English
 - ⊕ Simplified Chinese
 - ⊕ Spanish
 - ⊕ Russian

- ⊕ French
- ⊕ Italian
- ⊕ German
- ⊕ Norwegian
- Freeze Response
 - ⊕ Cine
 - ⊕ Comment
 - ⊕ Measure
 - ⊕ Body Mark
 - ⊕ Arrow
 - ⊕ Distance
- Trackball Sensitivity
 - ⊕ 1, 2, 3, 4, 5, 6
- Date/Time Setting
- Date Format
 - ⊕ MM/DD/YYYY, DD/MM/YYYY, YYYY/MM/DD
- Monitor Type
 - ⊕ TV-NTSC, TV-PAL, VGA
- Caps Lock
 - ⊕ ON, OFF
- Clear Annot On Unfreeze
 - ⊕ ON, OFF
- Guide Line Type
 - ⊕ Gun, Needle
- Focal Auto
 - ⊕ On, Off
- Display
 - Color of ROI
 - ⊕ Green, Yellow, Orange, Cyan
 - Display Format
 - ⊕ V1/2, V1/3, V2/3, H1/2, H1/4, O1/4
 - Screen Saver
 - ⊕ ON, OFF

- Screen Saver Delay (min)
 - ⊕ 0-99 minutes
- Annot Font Size
 - ⊕ Large, Medium or Small
- Storage
 - Clip Format
 - ⊕ AVI, WMV
 - Still Format
 - ⊕ BMP, JPG, TIF
 - Store Frame Amount
 - ⊕ 100-1000 frames
 - Print & Store Region
 - ⊕ Full Screen, Image Area, Image & Result Area
 - Image Share Service
 - ⊕ ON, OFF
 - Store to UsbDisk
- Key Configuration
 - The saved image can be printed, saved to USB drive, DICOM storage server and DICOM printer server.
 - The saved cine can be saved to USB drive and DICOM storage server.
 - P1 shortcut key: can be set to save image, cine or volume data or set to print saved image, send saved image or data to USB drive, DICOM storage server or DICOM printer server.
- Peripheral
 - Local Network
 - DHCP or Static IP
 - Enable DHCP, DICOM AE title and Mac address should be provided.
 - Enable Static IP, IP address, netmask, default gateway, DICOM AE title and Mac

- address should be provided.
- Advance
 - ⊕ Speed (10M, 100M, 1000M UNKNOW!)
 - ⊕ Duplex (Semi Duplex, Full Duplex, Unknow!)
- Ping IP Address function
- Printer
 - Printer Driver
 - Video Invert
 - Default Printer
 - Net Printer IP
- Comment
 - Classified by application parts
 - Support user-defined annotations
 - Support deleting user-defined annotations
 - Sort annotations
- Body Mark
 - Classified by application parts
 - Up to 135 types
- Measure
 - General
 - BSA
 - ⊕ Eastern
 - ⊕ Western
 - Cross Cursor Size
 - ⊕ Small, Medium, Large
 - Measure Line Size
 - ⊕ Small, Medium, Large
 - Dist Dash Line Disp.
 - ⊕ On, Off
 - Vel. Cross Line Disp
 - ⊕ On, Off
 - Ellip. Cross Line Disp

- ⊕ On, Off
- Line ID Disp.
- ⊕ On, Off
- Keep Result Window
- ⊕ On, Off
- Freeze on Measures
- ⊕ Yes, No
- Shortcut keys
- ⊕ Support performing obstetrics and gynecology measurements with shortcut keys. 10 shortcut keys can be set for each application.
- Formula
 - Fetal weight estimation
 - EFW GP/SD
 - Age by EFW
 - CUA
 - BPD
 - HC
 - FL
 - Cereb
 - GS
 - RAD
 - TIB
 - BOD
 - OFD
 - AC
 - HL
 - CRL
 - Clav.
 - Ulna
 - FIB
- Report
 - Logo

- Operation Logo (connect the USB drive written with report logo (170*60 pixels, bmp format) to the ultrasound system and click the left or right arrow key to select the desired logo)
- Title/Font
 - Title 1/2/3 (Font size can be set to 14, 16, 18, 20, 22, 24 and 26.)
 - Context (Font size can be set to 12, 13, 14, 15, 16)
 - Title (Font size can be set to 14, 16, 18, 20, 22, 24 and 26.)
 - Subtitle (Font size can be set to 14, 16, 18, 20, 22, 24 and 26.)
 - Background Color
 - Text Color
- Display Items For Report
 - Patient Information
 - Exam Information
 - Measurements (optional)
 - Image (optional)
 - Graphs (optional)
 - Comments (optional)
 - Preview button
- DICOM
 - Image Storage
 - Add or delete server list
 - Remote Host Name, IP Address, DICOM AE Title, Port Number of DICOM storage server; Connect Timeout, Dimse Timeout, Repeat Count, Acse Timeout, Relate to CommitStore, Send After Every Image Stored and Send At End of Exam.
 - Commitment Storage
 - Add or delete server list
 - Remote Host Name, IP Address, DICOM

AE Title, Port Number of DICOM Commitment storage; Connect Timeout, Dimse Timeout, Repeat Count, Acse Timeout, Synchronously Receive N-EVENT-REPORT Message

■ Worklist

- Add or delete server list
- Remote Host Name, IP Address, DICOM AE Title, Port Number of DICOM worklist server; Connect Timeout, Dimse Timeout, Repeat Count, Acse Timeout, Max. Result.

■ MPPS

- Add or delete server list
- Remote Host Name, IP Address, DICOM AE Title, Port Number of DICOM MPPS server; Connect Timeout, Dimse Timeout, Repeat Count, Acse Timeout

■ Print

- Add or delete server list
- Remote Host Name, IP Address, DICOM AE Title, Port Number of DICOM print server; Connect Timeout, print parameter pre-settings
- Format, Priority, Medium, Orientation, Film Size, Film Destination, Magnification, Copies, Smoothing Type, Trim, Border, Empty Image, Color, Min Density, Max Density, Configure Information, Film Session Label

■ QR List

- Add or delete server list
- Remote Host Name, IP Address, DICOM AE Title, Port Number of DICOM QR List server; Connect Timeout, Dimse Timeout, Repeat Count, Acse Timeout, Max. Result.

- Load Default
 - Load
 - New
 - Export/Import
- System Information

System Parameters

- Frame rate: up to 870fps (Probe dependent)
- Grayscale Level: 256
- Transducer Elements: Max. 256

B Mode

- Gain: 1-255 adjustable
- Scan Depth: 42.9 cm (Probe dependent)
- Image Zoom, Showing magnification (0.8-10 times)
- TGC: 8 levels slider controls
- Image Inversion: Left and Right, Up and Down
- Panoramic imaging: achievable
- Compound imaging: off, 1, 2
- Focus: Up to 12, Focus span adjustable
- Frequency: 5 bands adjustable
- Chroma: 13 types selectable
- Adaptive image fusion: 15 types selectable
- μScan: 0, 2, 3, 7, 11
- Line Density: 3 levels adjustable (High/Medium/Low)
- Persistence: 0-95 selectable
- Biopsy Guide Function: On/Off
 - Biopsy lines angle adjustable
 - Biopsy lines offset adjustable
- Dynamic Range: 20-280 (Probe dependent)
- Grayscale Curve: 7 selectable
- Imaging width and position: adjustable

- Power: 1-100 adjustable, one step each
- Tissue acoustic characteristics:1400-1700
- LGC: gain compensation for left or right part of image
- Trapezoid Imaging: On/Off (Linear probe)
- B steer mode (linear probe)
- M-tuning function

Color Doppler/TDI Mode

- Gain: 0-255
- Frame Rate: up to 50 frames/sec (probe dependent)
- Size and position of color ROI: adjustable
- Auto Focus (focus number:1)
- Inversion: Up/Down, Left/Right
- Flow Invert: On/Off
- Frequency Range: 5 steps, adjustable
- Wall Filter: 25-750Hz, adjustable (Probe dependent)
- PRF: 0.5-12kHz (Probe dependent)
- Line Density: 4 levels (low/medium/high/highest)
- Color/Direction energy: 10 levels for color Doppler mode; 4 levels for TDI mode
- Color baseline adjustment: ± 15 steps
- Persistence: 0-80 (Probe dependent)
- B Reject: 0-255
- Linear steer angle: 0, ± 16 , ± 20 adjustable
- M-tuning function

M Mode

- Steer M mode: 3 sample lines, Display frame rate
- Video Inversion (On/Off)

- Chroma: 5 types
- Display Format: H1/2, H1/4, V1/3, V1/2, V2/3, O1/4
- Scan Speed: 6 levels adjustable
- M Process: Switch between average and peak values
- Power: 30-100 adjustable

Spectral Doppler Mode

- Doppler Methods
 - PW (pulsed wave) Doppler
 - CW (continuous wave) Doppler
- Triplex: On/Off
- Sample Volume and Position for PW Doppler: 0.7-21mm adjustable
- Video Inversion: On/Off
- Spectrum Inversion: Achievable
- Angle Correction: 0°, 60°, -60°
- θ Angle Correction: correction range: 0 -72°
- Spectral Real-time Trace: Achievable
- Baseline Shift:17 steps selectable
- Frequency Range: 5 steps
- Wall Filter: 25-750 adjustable
- PRF: 1-16kHz (PW) (Probe dependent)
- PRF: 1-48KHz (CW) (Probe dependent)
- Max. Velocity Range
 - 0.0004-18 m/s (PW) (Probe dependent)
 - 0.0013-63 m/s (CW) (Probe dependent)
- Scan Speed: 4 levels adjustable
- Doppler Chroma: 5 kinds selectable
- One-key Auto Optimization
 - Auto Adjusting Baseline
 - Auto Adjusting PRF

- Auto Correcting Angle
- Dynamic Range: 10 types selectable
- Display Format: H1/2, H1/4, V1/3, V1/2, V2/3, O1/4
- Steer Angle: 0, ± 16 , ± 20 , 5 levels adjustable

3D Mode

- Display Mode
 - Dual-split display
 - Quad-split display
 - 3D Full Display
- Slice position: adjustable
- VOI feature: ON/OFF
- ROI size adjustable
- ROI sample line curvature: adjustable
- Restore: Volume/All/Curve, selectable
- Scan Mode: Lin and Sec, selectable
- Z Scale: 0.3-2.0 adjustable (range reduced for few slices)
- Z Angle: 10°-170° adjustable
- View angle: 0°, 90°, 180°, 270°
- Multi-Slice: A, B, C or 3D
- Render Mode: Surface, Skeleton, X-ray
- Crop plane: ON/OFF
- X Rotation
- Y Rotation
- Z Rotation
- Horizontal Movement: Left/Right
- Vertical Movement: Up/Down
- Rotate image freely by moving trackball
- Zoom In/Out: adjustable
- Transparency: 0-100 adjustable
- Contrast: 0-100 adjustable

- Brightness: 0-100 adjustable
- Smoothness: 0-30 adjustable
- Methods for cropping reviews
 - By trace (crop inner or outer image)
 - By box (crop inner or outer image)
 - By eraser (big or small eraser)
- Undo cut: undo last cut or restore the review
- Redo cut
- Auto rotate: 0°, 45°, 90°, 180°, 270°, 360° selectable
- Orientation: Top, Bottom, Left, Right, Front or Back, selectable
- 3D Chroma: 0-12 adjustable
- B Chroma: 1-13 adjustable
- Reference view: A, B or C
- Clip plane: AB, AC, BC or ABC
- M-slice: A, B or C
- Slice number: 3*3, 4*4 or 5*5
- Slice Spacing: 0.5-2.0, adjustable
- Measurement: Distance, Area, Volume
- 3D Image/Volume Storage
- Print Function
- Customized preset
- Focus position: adjustable
- Sweep Angle: 5° - 75°
- Image Quality: High, Medium, Low

4D Mode

- Display Mode
 - Dual-split display
 - Quad-split display
 - 3D Full Display
 - 4D Full Display
- Sample Line Curvature: adjustable

- ROI Size and Position: adjustable
- Focus position: adjustable
- Scan angle : 5°-75°, adjustable
- Image Quality: High, Medium, Low
- Stability: On/Off
- Slice position: adjustable
- VOI feature: ON/OFF
- ROI size adjustable
- ROI sample line curvature: adjustable
- Restore: Volume/All/Curve, selectable
- View angle: 0°, 90°, 180°, 270°
- Multi-Slice: A, B, C or 3D
- Render Mode: Surface, Skeleton, X-ray
- X Rotation
- Y Rotation
- Z Rotation
- Horizontal Movement: Left/Right
- Vertical Movement: Up/Down
- Rotate image freely by moving trackball
- Zoom In/Out: adjustable
- Transparency: 0-100 adjustable
- Contrast: 0-100 adjustable
- Brightness: 0-100 adjustable
- Smoothness: 0-30 adjustable
- Methods for cropping reviews
 - By trace (crop inner or outer image)
 - By box (crop inner or outer image)
 - By eraser (big or small eraser)
- Undo cut: undo last cut or restore the review
- Redo cut
- Cine Playback: The range depends on the volume data.)
- Auto Rotate: 0°, 45°, 90°, 180°, 270° and 360°, selectable

- Orientation: Top, Bottom, Left, Right, Front or Back selectable
- 3D Chroma: 0-12 adjustable
- B Chroma: 1-13 adjustable
- Reference view: A, B or C
- Clip plane: AB, AC, BC or ABC
- M-slice: A, B or C
- Slice number: 3*3, 4*4 or 5*5
- Slice Spacing: 0.5-2.0, adjustable
- Measurement: Distance, Area, Volume
- 3D Image/Cine/Volume Storage
- Customized preset
- Print function
- VolPre feature

Physiological Signal Display

- ECG Pulse wave
- ECG Lead-three lead system
- ECG Gain: adjustable
- ECG Position: adjustable
- ECG Invert: On/Off
- R-Trigger: On/Off
 - Trigger Delay: adjustable
 - Frame Count: adjustable

Integrated Data Management System

- Hard Disk Memory Capacity: 500G
- USB ports: 3 (including 1 engineering port)

Image Storage and Playback

- Cine Loop: Up to 1000 frames in B mode
- Real time single or dual display, static or dynamic image storage

- The Stored Images can be viewed directly on PC.
- Clip Board Function (In frozen B mode)
- Doppler Cine Playback: Speed is adjustable; Sound can be played back.

DICOM Network Communication

- Storage: Directly transmits images with patient information to a DICOM file server
- Print: Images can be printed directly using a DICOM compatible printer
- DICOM Storage Commitment
- DICOM Worklist
- DICOM MPPS
- DICOM Q/R
- Medical digital images and communication DICOM 3.0 interface

Preset Function

- Users can customize the presets based on different probe and diagnostic part to optimize imaging parameters and adjustment combination.
- Preset sequence adjustable
- Import/export the preset

Patient Data Management

- Patient Registration: Name, ID, Gender, Date of Birth, Height, Weight, LMP, EDD and GA.
- Patient data and reports are archived by patient exams
- Reports and images can be previewed.
- Image can be set to 1*1, 2*2, 3*3 for preview.
- Previewed file can be deleted, previewed, DICOM send or DICOM print.

- Data can be exported to USB drive or DVD in BMP, JPG, TIF, AVI or WMV format.

Annotation and Body Mark Setting

- Body Mark Icon: Classified by department
- Annotation can be selected and input in the library.

Physical Specification

- Size: 413mm*357 mm*391mm (L*W*D)
- Weight: Approx. 8.5kg
- Monitor: 15" Widescreen and High-Resolution Color LCD monitor, LED backlight, anti-flickering and vertically and horizontally rotatable
- Touch screen: 13.3"
- Probe connector: 2, interchangeable

Safety Standard

Comply with IEC60601-1, Class I , BF type

Environmental Requirements

- Operation Environment
 - Temperature: +10°C to +40°C
 - Relative Humidity: 30% to 85% (non condensing)
 - Atmospheric pressure: 700 to 1060hPa
- Transportation and Storage Environment
 - Temperature: -20°C to +55°C
 - Relative humidity: 20%- 90% (non condensing)
 - Atmospheric Pressure: 700 to 1060hPa
- Power Supply
 - 110-240~, 2.7-1.2A

- Frequency: 50/60Hz

Optional Probes

- Phased Array Probes
 - 4P-A((Frequency: 1.0-5.4MHz, sweep angle:90 °)
 - 3P1 (Frequency: 1.0-6.0MHz, sweep angle: 88 °)
 - 5P2 (Frequency: 2.0-9.0MHz, sweep angle: 89 °)
 - 8P1 (Frequency: 4.0-12.0MHz, sweep angle: 85 °)
 - MPTEE (Frequency: 4.0-13.0MHz, sweep angle: 90 °)
 - MPTEE mini (Frequency: 4.0-13.0MHz, sweep angle: 90 °)
- Linear Probes
 - L742 (4.0-16.0 MHz)
 - L752 (4.0-16.0 MHz)
 - 10I2 (4.0-16.0 MHz)
 - L743 (4.0-16.0 MHz)
 - 10L1 (4.0-16.0 MHz)
 - LAP7 (3.0-15.0MHz)
- Convex Probes
 - 3C-A (1.0-7.0 MHz)
 - C353 (2.0-6.8 MHz)
 - C322 (2.0-6.8 MHz)
 - C354 (2.0-6.8 MHz)
 - C362 (2.4-5.5 MHz)
 - C542 (3.0-6.5MHz)
 - 6V3 (3.0-15.0 MHz)
 - EC9-5 (3.0-15.0 MHz)
 - C611 (4.0-13.0 MHz)
 - C613 (4.0-13.0 MHz)

- 6V1 (3.0-15.0MHz)
- Volume Probes
 - VC6-2 (2.0-6.8MHz)
- Pencil Probes
 - PWD2.0
 - CWD2.0
 - CWD5.0
- Biplane Probes
 - BCC9-5
 - BCL10-5

Measurements and Calculations

- Basic Measurements and Calculations

- B Mode

Measurement Collection	Measurement Item
Distance	Two-Point Measurement, Length Trace Measurement, Distance Ratio Measurement, %Stenosis Distance
Area	Trace Area Measurement, Ellipse Area Measurement, Area Ratio Measurement, %Stenosis Area
Volume	Three-distance Measurement, Distance + Ellipse Measurement
Angle	Three-Point Angle Measurement, Two-Line Angle Measurement

- M Mode

Measurement Collection	Measurement Item
Basic Measurements	Distance Measurement, Slope Measurement, %Stenosis Distance Measurement, Distance Ratio Measurement, Time Measurement, Heart Rate Measurement

- Color-Flow Mode

Measurement Collection	Measurement Item
Distance	Two-Point Measurement, Length Trace Measurement, Distance Ratio Measurement, %Stenosis Distance
Area	Trace Area Measurement, Ellipse Area Measurement, Area Ratio Measurement, %Stenosis Area
Volume	Three-distance Measurement, Distance + Ellipse Measurement
Angle	Three-Point Angle Measurement, Two-Line Angle Measurement
Color Flow	Doppler Area Measurement, Color Flow Measurement, Flow Velocity Measurement

- Spectral Doppler-Mode

Measurement Collection	Measurement Item
Basic Measurements	Velocity Measurement, Acceleration Measurement, Resistivity Index Measurement, Pulsatility Index Measurement, S/D Ratio Measurement, Auto Trace Measurement, Manual Trace Measurement, Time Measurement, Heart Rate Measurement

● Obstetrics Measurements and Calculations

■ B Mode

➤ Measurements

Measurement Collection	Measurement Item
Fetal Biometry	BPD, OFD, HC, AC, FL, HL, Vp, CM, Cereb, NF
Early Gest	CRL, GS, YS, BPD, FL, NT
Long Bones	HL, Clav., RAD, Ulna, TIB, FIB
Fetal Cranium	Cereb., CM, Va, Vp, BOD, IOD, HEM, c.s.p, NT, NF
AFI	Q1, Q2, Q3, Q4
Uterus	Length, Height, Width, Endo.Thickn., Cervix Length
Lt(Rt) Ovary	Length, Height, Width
Umbilical Vein	Diameter
Lt (Rt) Uterine A	Diameter

➤ Calculations

GA
EDD
EFW
AUA
CUA

■ M Mode

Measurement Collection	Measurement Item
Generic	Distance, Slope, Time, HR, Stenosis % Distance
FHR	FHR, AtrialFHR

■ Spectral Doppler Mode

Measurement Collection	Measurement Item
Ductus Art	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR
Ao	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD,HR
Lt(Rt) Carotid	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR
Lt(Rt) MCA	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR

Measurement Collection	Measurement Item
Umbilical Art.	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR
Umbilical Vein	TAmax
Ductus Ven	S, D, a, PVIV, HR
Lt(Rt) Uterine Art.	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR
FHR	FHR
SMA	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR
Celiac.A.	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR
IVC	S, D, S,a, PLI, PVIV
Lt (Rt) UmA	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, HR

● Gynecology Measurements and Calculations

■ B Mode

Measurement Collection	Measurement Item
Uterus	Length, Height, Width, Endo.Thickn., Cervix Length
Lt(Rt) Uterine A	Diameter
Lt(Rt) Ovary	Length, Height, Width
Lt(Rt) Follicle	Follicle
Fibroid	Fibroid

■ M Mode

Measurement Collection	Measurement Item
Lt(Rt) Ovarian Art.	Vessel Diameter, Stenosis Diameter, Time, HR
Lt(Rt) Uterine Art.	Vessel Diameter, Stenosis Diameter, Time, HR
FHR	FHR, Atrial FHR

■ Spectral Doppler Mode

Measurement Collection	Measurement Item
Lt(Rt) Ovarian Art.	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, Time, HR
Lt(Rt) Uterine Art.	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED,

Measurement Collection	Measurement Item
	RI, SD, Time, HR
Vessel	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, Time, HR
FHR	FHR

● Cardiology Measurements and Calculations

■ B Mode

Measurement Category	Measurement Collection	Measurement Item
Dimensions	None	IVSd, LVIDd, LVPWd, IVSs, LVIDs, LVPWs, RVOT, AO, LA, ACS, LVOT, RVAWd, RVIDd, EPSS, MCS, MV Diam, AVA, MVA, MPA, PV Diam, TV Diam
Volume	Simp(LV)	A4Cd, A4Cs, A2Cd, A2Cs
	A-L(LV)	LVd, LVs
	Simp(LA)	A4Cs, A2Cs
	Simp(RA)	A4Cs
	Teichlozs(LV)	IVSd, LVIDd, LVPWd, IVSs, LVIDs, LVPWs
LV Mass	LV Mass (Cube)	IVSd, LVIDd, LVPWd
	LV Mass (A-L)	LVA _d Sax Epi, LVA _d Sax Endo, LVL _d Apical
	LV Mass (T-E)	LVA _d Sax Epi, LVA _d Sax Endo, a, d

■ Color-Flow Mode

Measurement Category	Measurement Collection	Measurement Item
Dimensions	None	IVSd, LVIDd, LVPWd, IVSs, LVIDs, LVPWs, RVOT, AO, LA, ACS, LVOT, RVAWd, RVIDd, EPSS, MCS, MV Diam, AVA, MVA, MPA, PV Diam, TV Diam
Volume	Simp(LV)	A4Cd, A4Cs, A2Cd, A2Cs
	A-L(LV)	LVd, LVs
	Simp(LA)	A4Cs, A2Cs
	Simp(RA)	A4Cs
	Teichlozs(LV)	IVSd, LVIDd, LVPWd, IVSs, LVIDs, LVPWs
LV Mass	LV Mass (Cube)	IVSd, LVIDd, LVPWd

Measurement Category	Measurement Collection	Measurement Item
	LV Mass (A-L)	LVAd Sax Epi, LVAd Sax Endo, LVLd Apical
	LV Mass (T-E)	LVAd Sax Epi, LVAd Sax Endo, a, d
PISA	PISA-MR	MR Rad, MR VTI (gray in the Doppler mode)
	PISA-AR	AR Rad, AR VTI (gray in the Doppler mode)
	PISA-TR	TR Rad, TR VTI (gray in the Doppler mode)
	PISA-MS	MS Rad, MS VTI (gray in the Doppler mode)

■ M Mode

Measurement Category	Measurement Collection	Measurement Item
Dimensions	None	IVSd, LVIDd, LVPWd, IVSs, LVIDs, LVPWs, RVOT, AO, LA, ACS, LVOT, RVAMd, RVIDd, EPSS, MCS
Time/Slope	/	LVPEP, LVET, RVPEP, RVET, MV DE, MV E-F Slope, MV E Amp, MV A Amp
	LV TEI	MV C-O Dur, LVET
	RV TEI	TV C-O Dur, RVET
Volume	Teichlozs(LV)	IVSd, LVIDd, LVPWd, IVSs, LVIDs, LVPWs
	Cube(LV)	IVSd, LVIDd, LVPWd, IVSs, LVIDs, LVPWs
LV Mass	LV Mass	IVSd, LVIDd, LVPWd

■ Spectral Doppler Mode

Measurement Category	Measurement Collection	Measurement Item
MV	None	MV E Vel, MV A Vel, MV E Dur, MV A Dur, MVA(PHT), MV DecT, MR Vmax, MR VTI
	MVA(VTI)	LVOT, LVOT VTI, MV VTI

MV	LV TEI	MV C-O Dur, LVET
	PISA-MR	MR Rad, MR VTI
	PISA-MS	MS Rad, MS VTI
AV	None	AV Vmax, AV VTI, AR Vmax, AR VTI, LVET, LVPEP, IVCT, IVRT, AR DecT, AR PHT, HR-AV
	AVA(VTI)	LVOT, LVOT VTI, AV VTI
	AVA(Vmax)	LVOT, LVOT Vmax, AV Vmax
	PISA-AR	AR Radius, AR VTI
TV	None	TV Vmax, TV VTI
	E/A	TV E Vel, TV A Vel
	RV TEI	TV C-O Dur , RVET
	RVSP	TR Vmax, RAP
	PISA-TR	TR Rad, TR VTI
PV	None	PR Vmax, PR VTI, PV AccT, MPA Vmax, RPA Vmax, LPA Vmax, RVET, RVPEP, HR-PV
	PAEDP	PR Ved, RAP
	PVA(Vmax)	RVOT, RVOT Vmax, PV Vmax
	PVA(VTI)	RVOT, RVOT VTI, PV VTI
Pulm&Hep Vein	None	Pulm S Vel, Pulm S VTI, Pulm A Vel, Pulm D Vel, Pulm D VTI, Pulm A Dur, Pulm DecT, Hep S Vel, Hep D Vel, Hep A Ve, Hep A Dur
TDI	/	Sa Medial, Ea Medial, Aa Medial, Sa lateral, Ea lateral, Aa lateral

● Abdomen Measurements and Calculations

■ B Mode

Measurement Collection	Measurement Item
Liver	Length, Height, Width
Port Vein	Portal Vein Diameter, Flow Diameter
Gallbladder	Length, Height, Width, Wall, CBD
Pancreas	Duct., Head, Body, Tail
Spleen	Length, Height, Width

Measurement Collection	Measurement Item
Kidney	Length, Height, Width
Renal Aorta	Vessel Area, Stenosis Area, Vessel Diameter, Stenosis Diameter, Flow Diameter
Aorta	Vessel Area, Stenosis Area, Vessel Diameter, Stenosis Diameter, Flow Diameter
Bladder	Length, Height, Width

■ M Mode

Measurement Collection	Measurement Item
Renal Aorta	Vessel Diameter, Stenosis Diameter, Time, HR
Aorta	Vessel Diameter, Stenosis Diameter, Time, HR

■ Spectral Doppler Mode

Measurement Collection	Measurement Item
Renal Aorta	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, Time, HR
Aorta	Auto Trace, Manual Trace, PS, ED, RI, PI, PS, ED, RI, SD, Time, HR
Port Vein	Velocity, Time

● Vascular Measurement

■ B Mode

Measurement Category	Measurement Collection	Measurement Item
Carotid	Lt(Rt) Subclav A	Vessel Diam Vessel Area %Stenosis Diam %Stenosis Area IMT
	Lt(Rt) CCA	
	Lt(Rt) Bulb	
	Lt(Rt) ICA	
	Lt(Rt) ECA	
	Lt(Rt) Vertebral A	

Measurement Category	Measurement Collection	Measurement Item
UE Art	Lt(Rt) Innom A	Vessel Diam Vessel Area %Stenosis Diam %Stenosis Area IMT
	Lt(Rt) Subclav A	
	Lt(Rt) Axill A	
	Lt(Rt) Brach A	
	Lt(Rt) Rad A	
	Lt(Rt) Ulnar A	
	Lt(Rt) Sup Palm A	
	Lt(Rt) Deep Palm A	
UE Vein	Lt(Rt) Innom V	Vessel Diam Vessel Area %Stenosis Diam %Stenosis Area IMT
	Lt(Rt) Subclav V	
	Lt(Rt) Int Jugular V	
	Lt(Rt) Axill V	
	Lt(Rt) Ceph V	
	Lt(Rt) Basilic V	
	Lt(Rt) Brach V	
	Lt(Rt) Med Cub V	
	Lt(Rt) Rad V	
	Lt(Rt) Ulnar V	
LE Art	Lt(Rt) Com Iliac A	Vessel Diam Vessel Area %Stenosis Diam %Stenosis Area IMT
	Lt(Rt) Ext Iliac A	
	Lt(Rt) Int Iliac A	
	Lt(Rt) Com Fem A	
	Lt(Rt) SFA	

Measurement Category	Measurement Collection	Measurement Item
LE Art	Lt(Rt) PFA	Vessel Diam Vessel Area %Stenosis Diam %Stenosis Area IMT
	Lt(Rt) Popl A	
	Lt(Rt) Ant Tib A	
	Lt(Rt) Post Tib A	
	Lt(Rt) Peron A	
	Lt(Rt) Dors Ped A	
LE Vein	Lt(Rt) IVC	Vessel Diam Vessel Area %Stenosis Diam %Stenosis Area IMT
	Lt(Rt) Com Iliac V	
	Lt(Rt) Ext Iliac V	
	Lt(Rt) Int Iliac V	
	Lt(Rt) Com Fem V	
	Lt(Rt) SFV	
	Lt(Rt) PFV	
	Lt(Rt) Popl V	
	Lt(Rt) Ant Tib V	
	Lt(Rt) Post Tib V	
	Lt(Rt) Peron V	
	Lt(Rt) GSV Thigh	
	Lt(Rt) GSV Calf	
	Lt(Rt) LSV	

■ M Mode

Measurement Category	Measurement Collection	Measurement Item
Carotid	Lt(Rt) Subclav A	Vessel Diam %Stenosis Diam Time HR
	Lt(Rt) CCA	
	Lt(Rt) Bulb	
	Lt(Rt) ICA	
	Lt(Rt) ECA	
	Lt(Rt) Vertebral A	
UE Art	Lt(Rt) Innom A	Vessel Diam %Stenosis Diam Time HR
	Lt(Rt) Subclav A	
	Lt(Rt) Axill A	
	Lt(Rt) Brach A	
	Lt(Rt) Rad A	
	Lt(Rt) Ulnar A	
	Lt(Rt) Sup Palm A	
	Lt(Rt) Deep Palm A	
UE Vein	Lt(Rt) Innom V	Vessel Diam %Stenosis Diam Time HR
	Lt(Rt) Subclav V	
	Lt(Rt) Int Jugular V	
	Lt(Rt) Axill V	
	Lt(Rt) Ceph V	
	Lt(Rt) Basilic V	
	Lt(Rt) Brach V	
	Lt(Rt) Med Cub V	

Measurement Category	Measurement Collection	Measurement Item
UE Vein	Lt(Rt) Rad V	Vessel Diam
	Lt(Rt) Ulnar V	%Stenosis Diam Time HR
LE Art	Lt(Rt) Com Iliac A	Vessel Diam %Stenosis Diam Time HR
	Lt(Rt) Ext Iliac A	
	Lt(Rt) Int Iliac A	
	Lt(Rt) Com Fem A	
	Lt(Rt) SFA	
	Lt(Rt) PFA	
	Lt(Rt) Popl A	
	Lt(Rt) Ant Tib A	
	Lt(Rt) Post Tib A	
	Lt(Rt) Peron A	
	Lt(Rt) Dors Ped A	
LE Vein	Lt(Rt) IVC	Vessel Diam %Stenosis Diam Time HR
	Lt(Rt) Com Iliac V	
	Lt(Rt) Ext Iliac V	
	Lt(Rt) Int Iliac V	
	Lt(Rt) Com Fem V	
	Lt(Rt) SFV	
	Lt(Rt) PFV	
	Lt(Rt) Popl V	
	Lt(Rt) Ant Tib V	

Measurement Category	Measurement Collection	Measurement Item
LE Vein	Lt(Rt) Post Tib V	Vessel Diam %Stenosis Diam Time HR
	Lt(Rt) Peron V	
	Lt(Rt) GSV Thigh	
	Lt(Rt) GSV Calf	
	Lt(Rt) LSV	

■ Spectral Doppler Mode

Measurement Category	Measurement Collection	Measurement Item
Carotid	Lt(Rt) Subclav A	Auto Trace Manual Trace PS ED RI PI PS, ED, RI, SD HR
	Lt(Rt) CCA	
	Lt(Rt) Bulb	
	Lt(Rt) ICA	
	Lt(Rt) ECA	
	Lt(Rt) Vertebral A	
UE Art	Lt(Rt) Innom A	Auto Trace Manual Trace PS ED RI PI PS, ED, RI, SD HR
	Lt(Rt) Subclav A	
	Lt(Rt) Axill A	
	Lt(Rt) Brach A	
	Lt(Rt) Rad A	
	Lt(Rt) Ulnar A	
	Lt(Rt) Sup Palm A	
	Lt(Rt) Deep Palm A	

Measurement Category	Measurement Collection	Measurement Item
UE Vein	Lt(Rt) Innom V	Auto Trace Manual Trace PS ED RI PI PS, ED, RI, SD HR
	Lt(Rt) Subclav V	
	Lt(Rt) Int Jugular V	
	Lt(Rt) Axill V	
	Lt(Rt) Ceph V	
	Lt(Rt) Basilic V	
	Lt(Rt) Brach V	
	Lt(Rt) Med Cub V	
	Lt(Rt) Rad V	
	Lt(Rt) Ulnar V	
LE Art	Lt(Rt) Com Iliac A	Auto Trace Manual Trace PS ED RI PI PS, ED, RI, SD HR
	Lt(Rt) Ext Iliac A	
	Lt(Rt) Int Iliac A	
	Lt(Rt) Com Fem A	
	Lt(Rt) SFA	
	Lt(Rt) PFA	
	Lt(Rt) Popl A	
	Lt(Rt) Ant Tib A	
	Lt(Rt) Post Tib A	
	Lt(Rt) Peron A	
Lt(Rt) Dors Ped A		

Measurement Category	Measurement Collection	Measurement Item
LE Vein	Lt(Rt) IVC	Auto Trace Manual Trace PS ED RI PI PS, ED, RI, SD HR
	Lt(Rt) Com Iliac V	
	Lt(Rt) Ext Iliac V	
	Lt(Rt) Int Iliac V	
	Lt(Rt) Com Fem V	
	Lt(Rt) SFV	
	Lt(Rt) PFV	
	Lt(Rt) Popl V	
	Lt(Rt) Ant Tib V	
	Lt(Rt) Post Tib V	
	Lt(Rt) Peron V	
	Lt(Rt) GSV Thigh	
	Lt(Rt) GSV Calf	
Lt(Rt) LSV		

● **Urology Measurement**

■ B Mode

Measurement Collection	Measurement Item
Lt(Rt) Kidney	Length, Height, Width, Renal Cortex
Lt(Rt) Renal A	Vessel Diam, Vessel Area, %Stenosis Diam, %Stenosis Area
Lt(Rt) Vessel	Vessel Diam, Vessel Area, %Stenosis Diam, %Stenosis Area
Lt(Rt) Suprarenal	Length, Height, Width
Bladder	Length, Height, Width

Measurement Collection	Measurement Item
Prostate	Length, Height, Width
Lt(Rt) Testicle	Length, Height, Width

- M Mode

Measurement Collection	Measurement Item
Lt(Rt) Renal A	Vessel Diam, %Stenosis Diam, Time, HR
Lt(Rt) Vessel	Vessel Diam, %Stenosis Diam, Time, HR

- Spectral Doppler Mode

Measurement Collection	Measurement Item
Lt(Rt) Renal A	Auto Trace; Manual Trace; PS; ED; RI; PI; PS, ED, RI, SD; HR
Lt(Rt) Vessel	Auto Trace; Manual Trace; PS; ED; RI; PI; PS, ED, RI, SD; HR

- Small Part Measurement

- Thyroid

- B Mode

Measurement Collection	Measurement Item
Lt(Rt) Thyroid	Length, Height, Width
Lt(Rt) Sup. ParThyroid	Length, Height, Width
Lt(Rt) Inf. ParThyroid	Length, Height, Width
Thyroid Ist.	Length, Height, Width
Lt(Rt) STA	Vessel Diam
Lt(Rt) ITA	

➤ M Mode

Measurement Collection	Measurement Item
Lt(Rt) STA	Vessel Diam, Time, HR
Lt(Rt) ITA	Vessel Diam, Time, HR

➤ Spectral Doppler Mode

Measurement Collection	Measurement Item
Lt(Rt) STA	Auto Trace; Manual Trace; PS; ED; RI; PI; PS, ED, RI, SD; Time; HR
Lt(Rt) ITA	Auto Trace; Manual Trace; PS; ED; RI; PI; PS, ED, RI, SD; Time; HR

■ Testicle

➤ B Mode

Measurement Collection	Measurement Item
Lt(Rt) Testicle	Length, Height, Width

➤ M Mode

Measurement Collection	Measurement Item
Lt(Rt) Vessel	Vessel Diam, Time, HR

➤ Spectral Doppler Mode

Measurement Collection	Measurement Item
Lt(Rt) Vessel	Auto Trace; Manual Trace; PS; ED; RI; PI; PS, ED, RI, SD; Time; HR

- Breast

- B Mode

Measurement Collection	Measurement Item
Lt(Rt) Lesion1-5	Nipple-Les. Dist, Skin-Les. Dist, Length, Width, Height

- M Mode

Measurement Collection	Measurement Item
Vessel	Vessel Diam, Stenosis Diam, Time, HR

- Spectral Doppler Mode

Measurement Collection	Measurement Item
Vessel	Auto Trace; Manual Trace; PS; ED; RI; PI; PS, ED, RI, SD; Time; HR

- Pediatrics Measurement

- B Mode

Measurement Collection	Measurement Item
Hip	Hip Angle, Hip d-D Ratio

- Calculation Reports

- Reports for each departments

Classification	Descriptions
Vascular report	/
Obstetrics report	Measurement report Fetal growth curve report Fetal anatomical survey report Fetus compare (quadruplet) Comments

Classification	Descriptions
Gynecology report	/
Small part report	/
Cardiac report	/
Urology report	/
Abdomen report	/
Pediatrics report	/

- Measurement values in the report are editable.
- Image: selectable
- Information such as hospital name can be preset.
- Report can be output as HTML format.

NOTE:

- The specifications of this system may change without any prior notification.
- Some products or features may not be available in some countries.
- Please contact your local SonoScape sales representative for more information.

Service Information:

Address: 4/F, 5/F, 8/F, 9/F & 10/F, Yizhe Building, Yuquan Road, Nanshan, Shenzhen, 518051, Guangdong, China

Zip code: 518051

Tel: +86-755-26722890

Fax: +86-755-26722850

Email: service@sonoscape.net